

Sizing the Seattle Splint

1. Measure for an accurate size. Allow for a hand's width between back of knee and top of splint. If fitting the Seattle Splint to attach to the thigh extension, measure from the base of the heel to the crease behind the knee (popliteal space) with the leg extended.
2. Toes should not extend over the splint.
3. The splint must allow full knee flexion.
4. Do not over tighten the straps. There should be space for a finger to slide between the skin or dressing and the strap pad.

Sizing and Specifications

Size	Calf Length	Height	Unisex Shoe Size
100 / 2100 Extra Large	20" – 23"	6'1½" and taller	13 and larger
101 / 2101 Large	18" – 20"	5'10" – 6'1"	10 – 12.5
102 / 2102 Medium	16" – 18"	5'5" – 5'9½"	7 – 9.5
103 / 2103 Small	14" – 16"	5' – 5'4½"	5 – 7
104 Youth	12½" – 14"	3'5" – 5'	3 – 5
105 Child	9" – 12"		
106 Toddler	5" – 8"		

Disclaimer: The Seattle Splint and components are designed to be sized and customized by a qualified practitioner. PFS Med, Inc does not express or imply responsibility for injury or impairment in skin integrity resulting from splints improperly fit or inadequately supervised.

